

Monoclonal Antibodies

Murine Antibody Clone 39-15

Catalogue #	MA-5010
Specificity	Ventricular/slow Human Skeletal Myosin Light Chain-1 Amino Acids 34-44
Isotype	IgG1/k
Form	Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide
Purification	Protein G
Storage	For long term storage, aliquot and store at -20°C
Uses	ELISA, western blot, Immunoaffinity Purification
Notes	Matched monoclonal antibody (Cat. #MA-5030) available
Sizes Available	1.0 mg 5.0 mg 10.0 mg Bulk sizes and pricing available

For additional information or to request a quotation, please contact:
reagents@spectraldx.com

For in vitro research use only.