

# Monoclonal Antibodies

## Murine Antibody Clone 9FABP-60

Catalogue #	MA-6010
Specificity	Human Muscle Fatty Acid Binding Protein (15 kD)
Isotype	IgG1/k
Form	Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide
Purification	Protein G
Storage	For long term storage, aliquot and store at -20°C
Uses	ELISA, western blot
Notes	Matched Monoclonals (Cat. # MA-6030 and MA-6020) and Polyclonal (Cat. # PA-6010) antibodies available
Sizes Available	1.0 mg 5.0 mg 10.0 mg Bulk sizes and pricing available

For additional information or to request a quotation, please contact:  
[reagents@spectraldx.com](mailto:reagents@spectraldx.com)

**For in vitro research use only.**