

Monoclonal Antibodies

Murine Antibody Clone 1rFABP-299

Catalogue #	MA-6020
Specificity	Human Muscle Fatty Acid Binding Protein (15 kD)
Isotype	IgG1/k
Form	Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide
Purification	Protein G
Storage	For long term storage, aliquot and store at -20°C
Uses	ELISA
Notes	Matched Monoclonals (Cat. # MA-6010 and MA-6030) and Polyclonal (Cat. # PA-6010) antibodies available
Sizes Available	1.0 mg 5.0 mg 10.0 mg Bulk sizes and pricing available

For additional information or to request a quotation, please contact:

reagents@spectraldx.com

For in vitro research use only.

Spectral Medical Inc.

135 The West Mall, Unit 2 • Toronto, Ontario, Canada, M9C 1C2
Toll Free: 1.888.426.4264 • Tel: 1.416.626.3233 • Fax: 1.416.626.7383

ISO 9001/ISO-13485 Certified

spectraldx.com