

# Monoclonal Antibodies

SPECTRAL  
*Medical* INC

## Murine Antibody Clone BB-1F9

Catalogue #	MA-7010
Specificity	Human Glycogen Phosphorylase BB
Isotype	IgG1/k
Form	Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide
Purification	Protein G
Storage	For long term storage, aliquot and store at -20°C
Uses	ELISA, western blot, Immunoaffinity Purification
Reference	Clin. Chem., Vol. 45, No. 6, Supplement 1999
Sizes Available	1.0 mg 5.0 mg 10.0 mg  Bulk sizes and pricing available

For additional information or to request a quotation, please contact:  
[reagents@spectraldx.com](mailto:reagents@spectraldx.com)

**For in vitro research use only.**

Spectral Medical Inc.  
135 The West Mall, Unit 2 • Toronto, Ontario, Canada, M9C 1C2  
Toll Free: 1.888.426.4264 • Tel: 1.416.626.3233 • Fax: 1.416.626.7383

ISO 9001/ISO-13485 Certified

[spectraldx.com](http://spectraldx.com)