

Murine Antibody Clone 2CYTC-199

Catalogue # MA-8010

Specificity Human/Rat/Mouse native cytochrome C (14 kD)
Amino acids 60-63.
Does not react with Bovine/Porcine/Equine cytochrome C or apocytochrome C

Isotype IgG1/k

Form Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide

Purification Protein G

Storage For long term storage, aliquot and store at -20°C

Uses Affinity Purification, Immunoprecipitation, ELISA, Western Blot

Sizes Available 1.0 mg
5.0 mg
10.0 mg
Bulk sizes and pricing available

For additional information or to request a quotation, please contact:

reagents@spectraldx.com

For in vitro research use only.

Spectral Medical Inc.

135 The West Mall, Unit 2 • Toronto, Ontario, Canada, M9C 1C2
Toll Free: 1.888.426.4264 • Tel: 1.416.626.3233 • Fax: 1.416.626.7383

ISO 9001/ISO-13485 Certified

spectraldx.com