

Murine Antibody Clone 4CYTC-21

Catalogue # MA-8020

Specificity Cytochrome C (14 kD) denatured
Reacts with Human/Rat/Mouse/Bovine/Porcine/Equine
cytochrome C or apocytochrome C

Isotype IgG1/k

Form Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide

Purification Protein G

Storage For long term storage, aliquot and store at -20°C

Uses Western Blot

Sizes Available 1.0 mg
5.0 mg
10.0 mg
Bulk sizes and pricing available

For additional information or to request a quotation, please contact:
reagents@spectraldx.com

For in vitro research use only.