

Recombinant single chain cardiac Troponin I-C, second generation (ScTnI-C-2)

Catalogue #	RP-3700
Primary Structure	Human cTnI (aa28-110)-linker-cTnC
Source	E. coli
Purity	> 95% by SDS-PAGE
Molecular Weight	~30 kD by SDS PAGE
Concentration	Determined by Bradford Assay using IgG as the standard
Purification Method	Affinity purified
Storage Condition	-20 °C
Advantages	Contains no human or animal products
Uses	Calibrators and controls for immunoassays, western blot standards

Sizes Available	0.1 mg
	1.0 mg
Sizes Available	5.0 mg
	10.0 mg
	Bulk sizes and pricing available Quantities over 1.0 mg by special order.

Commercialization of RP-3700 requires licensing from Spectral Medical Inc.

For additional information or to request a quotation, please contact:

reagents@spectraldx.com

For research use only. Not for human or *in vivo* use.

Spectral Medical Inc.

135 The West Mall, Unit 2 • Toronto, Ontario, Canada, M9C 1C2
Toll Free: 1.888.426.4264 • Tel: 1.416.626.3233 • Fax: 1.416.626.7383

ISO 9001/ISO-13485 Certified

spectraldx.com