

## Streptavidin-NC (Streptavidin engineered to bind to nitrocellulose)

Catalogue #	RP-9000
Source	E. coli
Purity	> 90% by scanning SDS-PAGE
Molecular Weight	~ 24 kD for each subunit by SDS-PAGE
Purification Method	Affinity chromatography
Concentration	Determined by Bradford Assay using BSA as the standard
Storage Condition	-20 °C
Advantages	Contains no human or animal products
Uses	Calibrators and controls for immunoassays, western blot standards
Notes	This streptavidin molecule not only binds to nitrocellulose membrane readily but also preserves the full biotin binding ability.

Sizes Available	0.1 mg 1.0 mg 5.0 mg 10.0 mg  Bulk sizes and pricing available Quantities over 1.0 mg by special order.
-----------------	---

For additional information or to request a quotation, please contact:  
[reagents@spectraldx.com](mailto:reagents@spectraldx.com)

**For research use only. Not for human or *in vivo* use.**

Spectral Medical Inc.  
135 The West Mall, Unit 2 • Toronto, Ontario, Canada, M9C 1C2  
Toll Free: 1.888.426.4264 • Tel: 1.416.626.3233 • Fax: 1.416.626.7383

ISO 9001/ISO-13485 Certified

[spectraldx.com](http://spectraldx.com)