

Monoclonal Antibodies

Murine Antibody Clone 2CA-4

Catalogue #	MA-4010
Specificity	Human Carbonic Anhydrase III
Isotype	IgG1/k
Form	Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide
Purification	Protein G
Storage	For long term storage, aliquot and store at -20°C
Uses	ELISA, western blot
Sizes Available	1.0 mg 5.0 mg 10.0 mg Bulk sizes and pricing available

For additional information or to request a quotation, please contact:
reagents@spectraldx.com

For in vitro research use only.