

Monoclonal Antibodies



Murine Antibody Clone 2CA-2

Catalogue #	MA-4020
Specificity	Human Carbonic Anhydrase III
Isotype	IgG1/k
Form	Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide
Purification	Protein G
Storage	For long term storage, aliquot and store at -20°C
Uses	ELISA, western blot
Sizes Available	1.0 mg 5.0 mg 10.0 mg Bulk sizes and pricing available

For additional information or to request a quotation, please contact:
reagents@spectraldx.com

For in vitro research use only.

Spectral Medical Inc.
135 The West Mall, Unit 2 • Toronto, Ontario, Canada, M9C 1C2
Toll Free: 1.888.426.4264 • Tel: 1.416.626.3233 • Fax: 1.416.626.7383
ISO-13485 Certified

spectraldx.com