

Monoclonal Antibodies

Murine Antibody Clone 39-15

Catalogue # MA-5010

Specificity Ventricular/slow Human Skeletal Myosin Light Chain-1
Amino Acids 34-44

Isotype IgG1/k

Form Purified - Solution in PBS, pH 7.3 Does not contain Sodium Azide

Purification Protein G

Storage For long term storage, aliquot and store at -20°C

Uses ELISA, western blot, Immunoaffinity Purification

Notes Matched monoclonal antibody (Cat. #MA-5030) available

Sizes Available 1.0 mg
5.0 mg
10.0 mg
Bulk sizes and pricing available

For additional information or to request a quotation, please contact:

reagents@spectraldx.com

For in vitro research use only.