


Cardiac Biomarker Proteins (Reagents)

Single-Chain Cardiac Troponin I-C Second Generation (ScTnI-C-2)



Catalogue #	RP-3700
Primary Structure	Human cTnI (aa28-110)-linker-cTnC
Source	<i>E. coli</i>
Purity	> 95% by SDS-PAGE
Molecular Weight	~30 kD by SDS PAGE
Concentration	Determined by Bradford Assay using IgG as the standard
Purification Method	Affinity purified
Storage Condition	-20 °C
Advantages	Contains no human or animal products
Uses	Calibrators and controls for immunoassays, western blot standards
Sizes Available	0.1 mg 1.0 mg 5.0 mg 10.0 mg Bulk sizes and pricing available Quantities over 1.0 mg by special order.

Commercialization of RP-3700 requires licensing from Spectral Medical Inc.

For additional information or to request a quotation, please contact:

reagents@spectraldx.com

For research use only. Not for human or in vivo use.

Spectral Medical Inc.

135 The West Mall, Unit 2 • Toronto, Ontario, Canada, M9C 1C2

reagents@spectraldx.com or slahiri@spectraldx.com

Phone: 416-626-3233 Extn 2200.

SPECTRAL | **U**
Medical **DX**

spectraldx.com/cardiac-biomarker-proteins