


Cardiac Biomarker Proteins (Reagents)



Recombinant NT-ProBNP

Catalogue #	RP-8100
Primary Structure	Human NT-ProBNP (aa 1-76) plus His tag
Source	<i>E. coli</i>
Purity	> 90% by scanning SDS-PAGE
Molecular Weight	~ 9 kD for each subunit by SDS-PAGE
Purification Method	Affinity Chromatography
Concentration	Determined by Bradford Assay using BSA as the standard
Storage Condition	-20 °C
Advantages	Contains no human or animal products
Uses	Calibrators and controls for immunoassays, western blot standards
Sizes Available	0.1 mg 1.0 mg 5.0 mg 10.0 mg Bulk sizes and pricing available Quantities over 1.0 mg by special order.

For additional information or to request a quotation, please contact:

reagents@spectraldx.com

For research use only. Not for human or in vivo use.

Spectral Medical Inc.

135 The West Mall, Unit 2 • Toronto, Ontario, Canada, M9C 1C2

reagents@spectraldx.com or slahiri@spectraldx.com

Phone: 416-626-3233 Extn 2200.

SPECTRAL 
Medical

spectraldx.com/cardiac-biomarker-proteins